



POWER MOBILITY DEVICE (POWER WHEELCHAIR & SCOOTER) User Guide (2 Sided)

This guide applies to power wheelchairs and power scooters. Always follow your physician's instructions and the manufacturer's operating manual.

BEFORE EACH USE

- ✓ Ensure the battery is fully charged.
- ✓ Check that the joystick or tiller controls operate properly.
- ✓ Inspect tires, wheels, footrests, armrests, and seat for damage.
- ✓ Verify the freewheel (manual drive) lever is in the **Drive** position before operating.
- ✓ Confirm the pathway is clear of obstacles.

SAFE OPERATION

- Sit fully back in the seat with your feet on the footrests before turning the unit on.
- Operate the controls slowly until you are familiar with the device.
- Start, stop, and turn gradually.
- Always reduce speed before making turns.
- Drive only on firm, level surfaces whenever possible.
- Use caution when approaching doorways, ramps, thresholds, curbs, and elevators.

SAFETY PRECAUTIONS

- ⚠ Always wear your seat belt if provided.
- ⚠ Never exceed the manufacturer's maximum weight capacity.
- ⚠ Do not allow passengers to ride on the device.
- ⚠ Avoid steep inclines, loose gravel, wet grass, sand, snow, or ice.
- ⚠ Never drive through standing water.
- ⚠ Keep hands, clothing, and personal items away from moving wheels.
- ⚠ Turn the power OFF before transferring into or out of the chair or scooter.
- ⚠ Do not operate while under the influence of alcohol, drugs, or medications that impair alertness.

BATTERY CHARGING

- Charge the batteries every night or after each day of use.
- Plug the charger directly into a grounded electrical outlet.
- Charge only in a dry, well-ventilated area.
- Do not cover the charger while charging.
- Disconnect the charger before operating the device.
- Contact Prime Care Medical Equipment if batteries fail to hold a charge.

TRANSFERS

- ✓ Turn the power OFF before transferring.
- ✓ Engage wheel locks if applicable.
- ✓ Flip up or remove footrests and armrests only if instructed by your provider.
- ✓ Transfer slowly using prescribed transfer techniques.
- ✓ Ask for assistance if you feel weak or unsteady.

CLEANING & MAINTENANCE

- Wipe the chair or scooter with a damp cloth and mild detergent.
- Keep the joystick, controller, and battery connections dry.
- Check tire condition and inflation (pneumatic tires) regularly.
- Inspect for loose hardware or unusual noises.
- Do not attempt repairs. Contact Prime Care Medical Equipment for service.

IF YOUR DEVICE STOPS WORKING

- Verify the power is ON.
- Check the battery charge indicator.
- Ensure the charger is disconnected.
- Verify the freewheel lever is in the **Drive** position.
- If the problem continues, discontinue use and contact Prime Care Medical Equipment.

For any medical emergency, call 911 immediately.